Chromatography Chromatographie

Nexera—
UHPLC 3.0

UHPLC 3.0

UHPLC systems represent the latest development in liquid chromatography instrumentation. Together with the enhancements necessary to operate with sub 2 µm particle columns, the differentiation is based on versatility in handling various applications on the same instrument without changes.

LC-30A Nexera

The LC-30A Nexera UHPLC has been developed as an all-round LC for handling diverse liquid chromatography requirements, providing ultra fast, high-resolution analysis while reducing the user's workload.

The Nexera UHPLC system is based on the LC-20A prominence series platform. Nexera contains entirely

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new developments and many unique solution concepts incorporating Shimadzu's broad knowledge in development of analytical and testing equipment.

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Performance Range


Nexera provides a wide range of performance offering short analysis time and a maximum number of injections.

Depending on the main focus, the use of short columns is possible as well as use of standard columns for optimizing analytical time and resolution. Long columns (150 - 250 mm) filled with small particles (sub 2 µm) enable highest resolutions (see graph).

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Shim-Pack XR Columns

Containing 2.2 µm particles, Shim-Pack XR columns have been designed for fast analysis at pressures up to 100 MPa. They combine speed with moderate backpressure, high resolution and excellent reproducibility.

Available column types:

- Shim-Pack XR-ODS/ODS II
- Shim-Pack XR-Phenyl
- Shim-Pack XR-SIL
- Shim-Pack XR-ODS III (available as 2 mm ID with 150 and 200 mm length only)

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
LC-30AD Pump

The LC-30AD double micro piston pump covers conventional and ultra high pressure LC.

- small piston stroke (10 µL) for optimal gradient performance under virtually pulse-free conditions
- flow range 0.0001 mL/min - 3 mL/min with 130 MPa backpressure, from 3.0001 - 5 mL/min up to 80 MPa
- ideal for use with high sensitivity detection techniques: MS (MS/MS), fluorescence etc.
- controls the vacuum degasser DGU-20
- plunger rinsing kit as standard
- optional 2 channel solvent selection valve

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SIL-30AC Autosampler

In terms of carryover and reproducibility, the SIL-30AC autosampler exceeds the performance of its predecessor and strengthens its outstanding position in the market. It offers DIIMS injection, fixed-loop injection for ultra-fast separation with minimized peak width as well as overlapping injection, allowing a maximum reduction in cycle time.

- total volume injection type autosampler for operating pressures up to 130 MPa based on the DIIMS (direct injection with isolated metering system) principle
- injection volume range 0.1 to 50 µL as standard
- two-solution rinsing kit and three rinse solvents as standard
- needle internal/external multi-rinse function to ensure ultra-low carry over
- automatic pretreatment functions
- optional loop injection kit
- cooled 105 vial rack for 1.5 mL vials as standard

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Micro Reactor

Instead of the relatively large mixing chamber volume used for many years, a micro-reactor (MiRC) achieves much better results for standard applications with its small mixing volume of only 20 µL. A 180 µL mixer is available for mixing complex mobile phases.

- new generation of gradient mixing in HPLC and UHPLC based on micro reactor technology
- Standard Type 20 µL – low volume, high efficiency mixer for two solvents
- TFA Type 180 µL – high efficient mixer for two solvents, for cases where high-sensitivity absorbance detection is required for a mobile phase containing UV-absorbing substances such as trifluoroacetic acid (TFA)

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CTO-30A Column Oven

For high-temperature LC analysis up to 150 °C, the column oven offers efficient pre-heating for solvents to enable stable and reliable conditions. The post-column cooler reduces detector noise.

- column thermostat, block heating type
- temperature range RT +5 °C to 150 °C
- heat balancer for achieving uniform temperature over the column
- for max. two columns up to 150 mm length

- options: solvent pre-heater and post column cooler reducing detector noise
- high pressure FCV-32 AH switching valve (two positions, six port)

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
Nexera— UHPLC 3.0

To find your local Shimadzu contact please visit www.shimadzu.eu

Shimadzu Europa GmbH
Albert-Hahn-Str. 6-10 · D-47269 Duisburg

Tel.: +49 - (0)203 - 76 87-0
Fax: +49 - (0)203 - 76 66 25
Email: shimadzu@shimadzu.eu