Sum Parameters Summenparameter

TOC, TN_b
Laboratory Labor

TOC, TN_b Laboratory Labor

The TOC-V-family meets all the needs of modern laboratory TOC analysis – regarding performance, variety of systems and methods. Shimadzu's know-how and experience for decades have been integrated in the systems.

TOC-V_{WS} (Standalone)

The TOC-V_{WS} (standalone) and TOC-V_{WP} (PC version) systems carry out TOC analysis according to wet chemical oxidation. This is a powerful oxidation technique using combination of sodium persulphate and UV oxidation at 80 °C. These three features guarantee complete analysis of all dissolved carbon species.

- automatic reagent preparation mode
- large injection volumes (up to 20.4 mL)
- excellent reproducibilities in the ppb range, up to the ultra-trace level
- extremely low detection limit of 0.5 µg/L
- operation of the TOC-V_{WS} standalone version via the large integrated control panel display
- TOC-Control V software runs the TOC-V_{WP} PC version

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Summenparameter

ASI-V

Filtrated samples and suspensions can be analysed quickly and easily by connecting the ASI-V to the TOC-V series, with the exception of the TOC-V_E. Suspensions are homogenised within the device.

- high sample capacity of 93 or 68 vials, depending on sensitivity
- autosampler containing 24 sample vials (125 mL) for TOC trace analysis with the TOC-V_{WP/WS}

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TOC-V_{CSH/N}

TOC-V_{CSH/N} version with its internal software makes it possible to investigate the content of TOC in water samples quite easily. This can be done using the difference method ($TOC = TC - TIC$) as well as the addition method ($TOC = POC + NPOC$). It is also possible to analyse the TOC concentration in gas samples.

- catalytic method of oxidation at 680 °C with NDIR-Detection
- integrated sample preparation reduces workload
- device carries out acidification, dilution and sparging
- extended working range from 4 µg/L to 25,000 mg/L

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TOC-V_{CPH/N}

The TOC-V_{CPH/N} has the same properties as the TOC-V_{CSH/N}, however, it is controlled via the TOC-Control V software.

In this way, several processes in the areas of quality assurance and daily analyses, which are normally performed manually, can be assigned to the system. These include among other things control charts and the validation of methods according to DIN 38402, part 51.

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TOC-V_E

The basic model TOC-V_E allows TOC-analyses at a particularly favourable price.

- catalytic oxidation at 680 °C
- manual injection unit
- difference analytic ($TOC = TC - TIC$) and NPOC method

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TOC-V_{CSH} Online-Module (Option)

Using the TOC-V_{CSH} with this online module option, a sample stream can be measured continuously. The user has the entire specification of the TOC-V_{CSH} at his disposal. This module is designed for the online determination of ultrapure water up to surface water. The flexible stand-alone software meets all the requirements of a modern online system.

- automatic measurements of the USP/EP system suitability test
- supports offline measurements

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TNM-1

The TN unit enables a combined TOC/TN-analysis in a single measurement step.

- catalytic oxidation is generally set to 720 °C in this application
- Chemiluminescence detection for determining TN
- working range for the TN parameter covers up to 4,000 mg/L
- Detection limit 5 µg/L with TOC-V_{CPH/CSH}

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OCT-1

The OCT-1 offers a simple and cost-effective alternative to automation for a small number of samples. The sample vials can be used directly without the need for refilling. In addition, conventional homogenisation tools (such as magnetic stirrers) can be used.

The OCT-1 can be adapted to the entire TOC-V series with exception of the TOC-V_E.

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SSM-5000A

The system is complete with the solid sample module SSM-5000A.

- switch between aqueous and solid samples without further instrument conversions
- large sample size of up to 1 g
- separate TC determination (total carbon) and IC (inorganic carbon) enable the exact determination of the TOC contents
- high sensitivity enables its use in cleaning validation using the glass fibre swab method without the need for extraction

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TOC-Control V Software — CLASS Agent

The user-friendly TOC-Control V software operates TOC-V analyzers (with the exception of the TOC-V_E), visualizes the measuring sequence and records measured results according to GLP/GMP. All instruments in the TOC-V series – modern tools for straightforward and reliable TOC analysis in a broad range of application areas – are now enhanced with easy-to-use and attractive software features.

- clear layout of the user interface
- inserting additional samples during current analysis cycle
- recalculation of measured samples
- wide range of export and import capabilities
- change of log-in during instrument runtime
- supports 21 CFR Part 11 compliance requirements

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
CLASS Agent

- together with Shimadzu's CLASS Agent database, the TOC data can be acquired, archived and secured
- release of secured data through electronic signature

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TOC, TN_b Laboratory Labor

To find your local Shimadzu contact
please visit www.shimadzu.eu

Shimadzu Europa GmbH
Albert-Hahn-Str. 6-10 · D-47269 Duisburg

Tel.: +49 - (0)203 - 76 87- 0
Fax: +49 - (0)203 - 76 66 25
Email: shimadzu@shimadzu.eu